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## BONE MARROW ASPIRATION NEEDLE INSTRUCTIONS FOR USE:

ENG

**Indications For Use:** Intended for the purpose of harvesting bone marrow specimens.

**Contraindications:** For use only for aspiration of bone marrow as determined by a licensed physician. These sets should be used by a trained physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of bone marrow biopsy. Physician judgement is required when considering biopsy on patients with bleeding disorder, or receiving anti-coagulant medications.

**Possible Side Effects:** Pain, fever, bleeding, hematuria and infection.

**Caution:** Ethylen oxide sterilized. Non-pyrogenic. Check the package if package damaged or opened do not use the products. Check the expiry date and the gauge. Store in dry place, between 5 – 30 °C temperature, protect from sun light. Possible allergic reactions should be considered. It must be used only in hospitals by physicians who familiar with the device and procedures. For single Patient Use Only. Do not attempt to clean or resterilize this product. Do not use the device under MRI. After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

**Please double check for device function and possible damage for correct implantation and do not use if package is opened or damaged.**

### STERNAL ASPIRATION PROCEDURE:

1. Using the sterile technique, prepare the skin with antiseptic and drape.
2. Infiltrate the selected area with local anaesthetic.
3. Make a skin incision with a scalpel blade over the marked area.
4. Adjust the needle stop for correct needle length. The adjustable screw-stop allows for a minimum needle length of 15mm, when the screw-stop is removed, the needle length is 50mm.
5. Introduce the needle through the incision pointing toward the Sternum.
6. Using gentle, but firm pressure, advance the needle. Rotate the needle in an alternating clockwise - counterclockwise motion. The entrance into the marrow cavity is generally detected by a decreased resistance.
7. Remove the stylet by pulling the upper section of the handle straight out.
8. Attach a syringe with a luer taper to the hub of the biopsy needle with a firm, push, twist motion.
9. Apply suction by withdrawing syringe plunger. Remove the aspirated specimen collected from the syringe.


### ILIIAC CREST ASPIRATION PROCEDURE:

- Repeat steps 1 through 3 of the procedure and continue as follows:
4. The screw-stop can be removed for maximum needle length.
  5. Insert the needle into iliac crest in usual fashion.
  6. Using gentle, but firm pressure, advance the needle. Rotate the needle in an alternating clockwise - counterclockwise motion. The entrance into the marrow cavity is generally detected by a decreased resistance.
  7. Remove the stylet by pulling the upper section of the handle straight out.
  8. Attach a syringe with a luer taper to the hub of the biopsy needle with a firm, push, twist motion.
  9. Apply suction by withdrawing syringe plunger. Remove the aspirated specimen collected from the syringe.

**Note:** These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

Single Use

Storage limits

Expiry Date

STERILE EO

Ethylen oxide sterilized