

TUMOR MARKING NEEDLE TÜMÖR İŞARETLEME İĞNESİ

Sterile placement needle with preloaded 99.99% Pure Gold or PEEK® marker
Accurate target visualization in tumors/organs during external radiation therapy
Excellent visibility of markers under fluoroscopy, X-ray, CT, MRI and ultrasound
Knurled surface of the markers prevent migration



Handle #2

Distal ucuna altın ve PEEK® işaretleyicisi yüklenmiş olan steril iğne
Radyasyon tedavilerinde tümör veya organlarda kesin hedefi görüntüleme
Floroskopi, X-ray, CT, MRI ve ultrason altında mükemmel marker görünümü
Pürüzlü yüzeye sahip markerlar ile en iyi şekilde marker konumlandırılması sağlanır

PEEK MARKER

GOLD MARKER

Ref. No	Needle (G x cm)	PEEK® Marker (mm x mm)	Marker Shape	Ref. No	Size (G x cm)	Marker Size (mm x mm)	Marker Shape
Handle #2				Handle #2			
GTP1410#2	14 x 10	1,5 x 4	∞	GTM1815#2	18 x 10	1 x 4	⌋⌋⌋⌋⌋⌋
GTP1412#2	14 x 12	1,5 x 4	∞	GTM1815#2	18 x 15	1 x 4	⌋⌋⌋⌋⌋⌋
GTP1415#2	14 x 15	1,5 x 4	∞	GTM1820#2	18 x 20	1 x 4	⌋⌋⌋⌋⌋⌋
GTP1610#2	16 x 10	1,2 x 4	∞	GTM1825#2	18 x 25	1 x 4	⌋⌋⌋⌋⌋⌋
GTP1615#2	16 x 15	1,2 x 4	∞	GTM1915#2	19 x 15	0,8 x 4	⌋⌋⌋⌋⌋⌋
GTP1810#2	18 x 10	1 x 4	⌋⌋⌋⌋⌋⌋	GTM1920#2	19 x 20	0,8 x 4	⌋⌋⌋⌋⌋⌋
GTP1815#2	18 x 15	1 x 4	⌋⌋⌋⌋⌋⌋	GTM1925#2	19 x 25	0,8 x 4	⌋⌋⌋⌋⌋⌋
GTP1820#2	18 x 20	1 x 4	⌋⌋⌋⌋⌋⌋				
GTP1825#2	18 x 25	1 x 4	⌋⌋⌋⌋⌋⌋				
GTP1910#2	19 x 10	1 x 4	⌋⌋⌋⌋⌋⌋				
GTP1915#2	19 x 15	0,8 x 4	⌋⌋⌋⌋⌋⌋				
GTP1920#2	19 x 20	0,8 x 4	⌋⌋⌋⌋⌋⌋				
GTP1925#2	19 x 25	0,8 x 4	⌋⌋⌋⌋⌋⌋				

* Please indicate the type of handle (Handle #1 or Handle #2) when ordering.




Fig-1 Mammography capture of PEEK marker 1.2mmx5mm