











# DISPOSABLE NEEDLE GUIDES

COMPATIBLE USG PROBES		REFERENCE CODE	NEEDLE GUIDE
<b>HITACHI</b>	EUP-V33, EUP-V33W, EUP-CC331 EUP-CC531, EUP-V53W, EUP-V12	<b>GTK01</b>	
<b>BK</b>	8808	<b>GTK02</b>	
	8818	<b>GTK13</b>	
<b>TOSHIBA</b>	PVF-620ST, PVL-625RT PVM-620ST	<b>GTK03</b>	
	PVT-661VT, PVM-651VT, PVQ-641V, PVG-630V	<b>GTK04</b>	
	PVM-740RT, PVT-770-RT	<b>GTK07</b>	
	SSA-790A AplioXG	<b>GTK08</b>	
	PVT-781-VT, PVU-781VT PVT-781-VTE, PVU-781-VTE	<b>GTK130</b>	
<b>ESAOTE</b>	EC123, EC1123, SE3123	<b>GTK08</b>	
	E8-5R10	<b>GTK11</b>	
<b>ALOKA</b>	UST-670P, UST-675P, UST-676P,	<b>GTK06</b>	
	UST-984-5, UST-9112-5	<b>GTK100</b>	
	UST-9118, UST-9124	<b>GTK101</b>	
<b>SHIMADZU</b>	UBR10R-65U	<b>GTK05</b>	
<b>PHILIPS</b>	BP10-5EC C9-4EC	<b>GTK08-R</b>	
	C8-4V, C9-3V C9-4V	<b>GTK110</b>	

COMPATIBLE USG PROBES		REFERENCE CODE	NEEDLE GUIDE
<b>GE</b>	IC5-9RS, E72, E721, E8C	<b>GTK08</b>	
	E8CS, IC5-9-D IC5-9H	<b>GTK09</b>	
	6.5MHz	<b>GTK09-S</b>	
	BE9C, BE9CS	<b>GTK10</b>	
	E7C-RC	<b>GTK11-R</b>	
<b>ALPINION</b>	RIC-9W-RS, (Voluson S) RIC6-12-D, RIC5-9D, IC9-RS	<b>GTK140</b>	
	(VOLUSON S150)	<b>GTK141</b>	
<b>SIEMENS</b>	E3-10	<b>GTK08</b>	
	BE9-4 EC9-4 Acuson X-150	<b>GTK08</b>	
	EC7, EC10C5 6.5EC10, EC9-4 SONOLINE G20	<b>GTK11</b>	
<b>SAMSUNG MEDISON</b>	EC9-4W Acuson S3000 MC9-4	<b>GTK12</b>	
	Endo-V II Sonoline G50	<b>GTK120</b>	
	Endo V	<b>GTK121</b>	
	EC9-4 (Sonoline G40, Acuson X150 Acuson X300) EV9-4(Sonoline G40 Acuson X150, Acuson X300)	<b>GTK122</b>	
<b>MINDRAY</b>	EC4-9IS, EC4-9ES ER4-9ES	<b>GTK11</b>	
<b>MINDRAY</b>	V10-4B, V10-4 CB10-4, 6CV1s	<b>GTK11</b>	

\* Please contact us for your specific requirements or if you can not find the probe in the list.

# DISPOSABLE NEEDLE GUIDES



- ✓ Used with under ultrasound guided procedures for Transrectal prostate biopsies, Transvaginal biopsies and Transvaginal oocyte retrieval.
- ✓ The disposable designs reduce the risks of cross contamination and infection.
- ✓ All guides accept 16G -18G sizes instruments.
- ✓ Compatible with well-known endocavity ultrasound probes in the market. See table below.

**GEOTEK**  
Healthcare Products

İvedik OSB Ağaç-Metal San. Sit. 1436. Sokak No: 12 06370 Ostim - Ankara / TÜRKİYE  
Tel: +90 (312) 395 98 52 • Fax: +90 (312) 395 98 53  
e-mail: info@geotekmedical.com • www.geotekmedical.com

**GEOTEK**<sup>®</sup>  
Healthcare Products  
CE 1984